


Optima

**LAMIDEY
NOURY**
MEDICAL



The complete system for electrosurgery in one unit



Vessel sealing system

Safe and efficient sealing up to 7mm diameter



Thermocision

A surgery revolution ! The only system which can dissect, seal and cut




Saline cut bipolar

A specific current for each procedure : TURP / TURB / HYSTERO RESECTION



Monopolar / Bipolar currents

A real hemostatic cut with a white coagulation without charring and resection without bleeding



Argon system

Argon module for a complete solution

FRIENDLY USER

- Touch screen
- Keyboard touch screen
- 4 High Frequency outputs
- Programs library
- Saving personal programs
- Automatic bipolar
- Instructions videos saved for instruments setting
- Multi language

MAINTENANCE




- Upgradeable by software
- Communication link with computer

SAFETY

- PFB™ : Power Feed Back
- STAHF™ : Automatic HF stabilization
- Ready for use
- SENTM™ : smart safety management
- Power limit for neonatal surgery



Unit Range

RANGE	OPTIMA 4	OPTIMA 3	OPTIMA 2
			
POWER	370W	370W	200W
OUTPUTS	2 Monopolars 2 Bipolars	1 Monopolar 2 Bipolars	1 Monopolar 1 Bipolar
OPTIONS	Vessel sealing Saline cut bipolar	Vessel sealing Saline cut bipolar	

Reusable Vessel Sealing Forceps



Lamidey

Z.A. Les Godets 3 rue des Petits Ruisseaux
F 91370 Verrières le Buisson France
Tél. : +33 (0)1 69 20 69 69
Fax : +33 (0)1 60 13 97 47
info@lamidey-noury.fr
www.lamidey-noury.fr

French manufacturer of electrosurgical unit since 1947